



Hygie® Alu-Classic Commode Chair

Lightweight and hygienic. Simple and easy to clean.



No other chair is easier to CLEAN and MAINTAIN than the Hygie® Alu-Classic Commode Chair.

- ✓ 100% Aluminum
- ✓ Rustproof
- ✓ Lighter
- ✓ More manageable
- ✓ Easier to clean

Used with the Commode Hygienic Covers® with super-absorbent pad



- The Hygie® Alu-Classic Commode Chair is a commode with swivel wheels, independent brakes, and pivoting armrest. The Alu-Classic commode allows height adjustment for additional comfort for patients and a more comfortable posture for care staff: 15 in to 17 in or 18 in to 20 in depending on the model of the commode.
- The commode provides back support and has push type handles. Designed for positioning only, the Alu-Classic Commode Chair can replace the toilet by using Hygienic covers® with super-absorbent pad and can also be used over a standard toilet.
- Using Hygienic covers® with the Hygie® Alu-Classic Commode Chair creates the perfect solution to contain potentially harmful pathogens and further promote safe, reliable practices in healthcare. The unique super-absorbent pad included in each Hygienic cover® is designed to turn liquid into a gel within seconds.

www.hygie.com

Toll free : 1 (866) 588-2221

Fax : (450) 444-6222

Canada, USA



Hygie® Alu-Classic Commode Chair

Lightweight and hygienic. Simple and easy to clean.

	Models 15, 16 and 17 inches		Models 18, 19 and 20 inches	
Total width	57 cm	22.5 in	57 cm	22.5 in
Total length	76 cm	30.0 in	76 cm	30.0 in
Height adjustable seats	38 cm	15.0 in	47 cm	18.5 in
	40 cm	16.0 in	49.5 cm	19.5 in
	43 cm	17.0 in	52 cm	20.5 in
Low height handles	92.7 cm	36.5 in	101.6 cm	40 in
Seat dimension				
	Width	48.3 cm	19.0 in	48.3 cm
Depth	44.5 cm	17.5 in	44.5 cm	17.5 in
Maximum charge	400 lbs	181.4 kg	400 lbs	181.4 kg

**Up to
400 lbs**

● **Frame**

100% aluminum fabrication.

● **Surface**

Can be disinfected with common commercial solutions. It does not contain porous material.

● **Spare parts**

Hygie can supply all spare parts, whether they are under warranty or not.

● **High stability**

The Hygie® Alu-Classic Commode Chair has a wider wheelbase chassis than most commodes found on the market, which makes it safer and more stable during use.

● **Armrest swivel**

The arms can be lowered by pivoting on the axis of the chair, completely freeing the seat surface to facilitate transfers.

● **Capacity**

The maximum capacity of the Hygie Classic® Commode Chair is of 180 kg/400 lbs.

● **Swivel and adjustable wheels**

4-inch swivel wheels with total lock brakes that won't mark the floor. The adjustable wheels allow for 3 different heights. Floor-to-seat clearance for 18, 19 or 20 inches. The small commode chair option offers floor clearance of 15, 16 and 17 inches.

To order:

- CH-HYAL-SBCU-018 : Bedpan seat - Height 18"
- CH-HYAL-SBHY-018 : Bedpan seat with hole - Height 18"
- CH-HYAL-SBSF-018 : Bedpan seat without bottom - Height 18"
- CH-HYAL-SBCU-016 : Bedpan seat - Height 16"
- CH-HYAL-SBHY-016 : Bedpan seat with hole - Height 16"
- CH-HYAL-SBSF-016 : Bedpan seat without bottom - Height 16"

www.hygie.com

Toll free : 1 (866) 588-2221

Fax : (450) 444-6222

Canada, USA




 **hygie[®]**

ALU-CLASSIC
COMMODE CHAIR

**NEW
TECHNOLOGY**

“ **One** out of nine patients admitted to hospital each year develops hospital acquired infections”*



New technology in the prevention of infections

On the leading edge of the battle against nosocomial infections, Hygie is proud to introduce its brand new Alu-Classic® commode chair, made entirely of aluminium and free of porous surfaces. This Canadian-built product is specially designed for hospitals and long-term care facilities.

"This chair responds to infection control needs in health care facilities, as it contains no porous material. This helps facilitate cleaning and disinfection."

France - Infection Control Nurse

"The Alu-Classic commode chair is a must for hospitals and long-term care facilities. It's designed to facilitate access on both sides while providing patients with greater comfort and security."

Isabelle - Occupational therapist


"We are convinced that this new generation of commode chairs will significantly reduce patient falls. The rubberized wheels provide greater adherence and the four-wheel braking system ensures the chair remains completely motionless."

Claude - Nursing assistant



Technical specifications

1	Total width	57 cm	22.5 po
2	Total height	89 to 102 cm	35 in to 40 in
3	Total depth	68 cm	27 in
4 A	Floor to seat height (front)	38 to 43 cm	15 in to 17 in
4 B	Floor to seat height (front)	46 to 51 cm	18 in to 20 in
5	Floor to seat height (back)	37 to 49 cm	14.5 in to 19,5 in
6	Seatback	42 cm	16.5 in
7	Seatback angle	5°	
8	Seat width	49 cm	19 in
9	Seat depth	44 cm	17.5 in
10	Distance between handles	15 cm	5.75 in
11	Clearance height	34 to 47 cm	13.5 in to 18.5 in
12	Pivoting armrest height	20 cm	8 in
13	Wheel diameter	15 cm	6 in
14	Grip	13 cm	5 in
15	Side rail	51 cm	20 in
16	3 height adjustments	2.5 to 7.6 cm	1 in to 3 in
17	Floor to side rail height	20 cm	8 in
	Chair weight	9 kg	20,5 lb
	Maximum load	181 kg	400 lb



From 1990 to today

Over 80 options
to meet your
every need

HYGIE ALU-CLASSIC® CHAIR

- ✓ Aluminium construction
- ✓ Ultra-light, sturdy and stable
- ✓ Adjustable height
- ✓ Easy to clean
- ✓ Pivoting armrests
- ✓ Free of porous surfaces
- ✓ Rotates freely on itself
- ✓ Rust-free
- ✓ Floor to seat height of 15"



HYGIE CLASSIC® CHAIR

- ✓ Galvanized steel construction
- ✓ Sturdy and stable
- ✓ Adjustable height
- ✓ Easy to clean
- ✓ Pivoting armrests
- ✓ Free of porous surfaces
- ✓ Large size available
24" wide



WARRANTY: Chair frames are guaranteed for a period of two (2) years (Classic Chair) or five (5) years (Alu-Classic Chair) respectively from the date of purchase. Moving pieces: armrests, wheels, etc. are guaranteed for a period of one (1) year from the date of purchase. Accessories (boots (armrest receptacles), pails, seats, etc.) are guaranteed for a period of three (3) months from the date of purchase. In the case of parts or pieces with a manufacturer's defect, Hygie® will supply replacement parts of original, similar or equivalent quality.

Wheelbase

The wheelbase of the Hygie Alu-Classic commode chair's frame makes it highly stable and secure for the patient, reducing the risk of overturning.

Rubber wheels

The 4-inch rubberized wheels offer 360° pivoting action and are equipped with individual manual brakes to fully immobilize the chair. The rubber also improves floor grip without leaving scuff marks.

Adjustable height

The Alu-Classic commode chair allows height adjustment for additional comfort for patients and a more comfortable posture for care staff (see Technical Specs for models 4A and 4B on page 3).



360°



Pivoting armrests

The sturdy, secure armrests can be lowered by pressing a button. They also pivot along the seat axis, which frees that area and facilitates patient transfers.



Easy to clean

This chair was specially designed to minimize gaps and hard-to-reach spaces, ensuring quick and easy disinfection.



Construction

The lightness of the chair's frame is due to its aluminium construction. The slight tilt in the seatback provides better comfort to the patient.



Single-patient seats

Use of Hygienic Cover® recommended



1 AC-SIEG-HY21-000
Bedpan Seat with holes. Use with Hygie Hygienic Cover®

Cleaning bedpan seats and pails

- We recommend the use of detergents or enzyme products.
- Seats, pails and bedpans may be cleaned with a decontamination washer (maximum temperature of 90.5°C or 195°F).



2 CH-SIRO-OUSC-001
Round seat, open front without cover



3 AC-SIEG-HY21-SFO
Bottomless Bedpan Seat



4 CH-CUVE-HYAL-000
Pail with holes



5 CH-CUVE-HYAL-001
Pail with full bottom



6 AC-SIEG-CUVE-000
Bedpan Seat with bottom



Hygienic Cover®
SA-HYGI-CLBI-000



Removable accessory plate with Hygienic Cover® dispenser and a bracket to hold the disinfecting wipe receptacle. Install on the commode chair seatback.

Product code
CH-PLAQ-DIST-ALLU

The Hygie Hygienic Covers® for Commode Chairs are the perfect solution to effectively manage and contain body wastes at the point of care. They encourage healthcare providers to adopt optimal work practices, which reduce the spread of pathogens that causes epidemics, cross-contamination and nosocomial (hospital-acquired) infections.

The super-absorbent pad in each Hygienic Cover® is unique – it absorbs 700 ml of body fluids within 30 seconds and transforms them into a gel. Its use helps decrease the risk of splashes, spills and the spread of harmful germs, thus preventing nosocomial infections and cross-contamination.

Ease of maintenance and cleaning

- Designed to comply with new disinfection standards.
- To clean, use products with high detergent and disinfection properties such as a 5,000 ppm sodium hypochlorite solution following every use.
- Remove the bedpan prior to cleaning the frame.
- Begin the cleaning/disinfecting process, moving from least to most-soiled areas.
- We do not recommend the use of wheelchair washing machines. The mechanical action in manual cleaning is designed to enhance detergent efficiency by eliminating a high proportion of organic waste.
- Do not use or clean in the shower.
- Do not use a scouring pad.



Preventable statistics

10% of patients admitted into a short-term care facility within the healthcare system will contract an infection as a result of received care or treatments.⁽¹⁾

30% to 50% of hospital-acquired infections are preventable.⁽²⁾

250,000 Canadians admitted to hospital each year develop hospital-acquired infections and, of these patients, **8,000** die.⁽³⁾

Nosocomial infections cost the Quebec healthcare system approximately **\$180 million** annually.⁽⁴⁾

Nosocomial infections are the **fourth** leading cause of death in Canada.⁽⁵⁾

(1) <http://publications.msss.gouv.qc.ca>

(2) Canadian Committee on Antibiotic Resistance, 2008.

(3) Canadian Institutes of Health Research, 2012.

(4) <http://publications.msss.gouv.qc.ca>

(5) Statistics Canada. Ranking and number of deaths for the 10 leading causes, Canada, 2000 and 2009.





 9140 Leduc Blvd.,
Suite 230
Brossard (QC) J4Y 0L2
Canada

 www.hygie.com

 info@hygie.com

 1.866.588.2221

2915 Ogletown Road,
2658
Newark, DE 19713
U.S.A.

 1.866.588.2221

310, rue du Vallon,
imm l'Atlas - Les Bouillides
Valbonne, 06560
France

 008.600.008 #234







©2017, Hygie Canada inc. All rights reserved
Printed in Canada

Printed on recycled and recyclable paper
with vegetable ink



ALU-CLASSIC COMMUNE CHAIR

Place an order: po@hygie.com

COMMUNE OPTIONS	PRICE	HEIGHT	PRODUCT CODE
Hygie® Alu-Classic Hy21® - Bedpan seat with bottom Standard, 4" swivel casters with total lock Distance between armrests: 19.5"/495 mm Capacity: 181 kg/400 lb" Codes CHHYALSBCU015 to 17: 3 adjustments :15", 16" and 17" Codes CHHYALSBCU018 to 20: 3 adjustments :18", 19" and 20"		15" / 381 mm	CHHYALSBCU015
		16" / 406 mm	CHHYALSBCU016
		17" / 432 mm	CHHYALSBCU017
		18" / 457 mm	CHHYALSBCU018
		19" / 482 mm	CHHYALSBCU019
		20" / 508 mm	CHHYALSBCU020
Hygie® Alu-Classic - Bottomless bedpan seat Standard, 4" swivel casters with total lock Distance between armrests: 19.5"/495 mm Capacity: 181 kg/400 lb" Codes CHHYALSBSF015 to 17: 3 adjustments: 15", 16" and 17" Codes CHHYALSBSF018 to 20: 3 adjustments: 18", 19" and 20"		15" / 381 mm	CHHYALSBSF015
		16" / 406 mm	CHHYALSBSF016
		17" / 432 mm	CHHYALSBSF017
		18" / 457 mm	CHHYALSBSF018
		19" / 482 mm	CHHYALSBSF019
		20" / 508 mm	CHHYALSBSF020
Hygie® Alu-Classic - Hy21® Bedpan seat with holes Standard, 4" swivel casters with total lock Distance between armrests: 19.5"/495 mm Codes CHHYALS BHY015 to 17: 3 adjustments: 15", 16" and 17" Codes CHHYALS BHY018 to 20: 3 adjustments: 18", 19" and 20"		15" / 381 mm	CHHYALS BHY015
		16" / 406 mm	CHHYALS BHY016
		17" / 432 mm	CHHYALS BHY017
		18" / 457 mm	CHHYALS BHY018
		19" / 482 mm	CHHYALS BHY019
		20" / 508 mm	CHHYALS BHY020
Hygie® Alu-Classic with round open front seat without cover Standard, 4" swivel casters with total lock Distance between armrests: 19.5"/495 mm Capacity: 181 kg/400 lb Pail not included. Codes CHHYALROSC015 to 17: 3 adjustments: 15", 16" and 17" Codes CHHYALROSC018 to 20: 3 adjustments: 18", 19" and 20"		15" / 381 mm	CHHYALROSC015
		16" / 406 mm	CHHYALROSC016
		17" / 432 mm	CHHYALROSC017
		18" / 457 mm	CHHYALROSC018
		19" / 482 mm	CHHYALROSC019
		20" / 508 mm	CHHYALROSC020

5 year warranty on the chair frame. Any items that have been modified or stored, installed and used incorrectly are not covered by the warranty.



ACSIEGHY21000



ACSIEGHY21SF0



ACSIEGCUVE000



CHSIROOUSC001



CHCUVEHYAL000



CHCUVEHYAL001

ACCESSORIES

SEATS	PRICE	UNITS PER CASE	PRODUCT CODE
Hy21° bedpan seat with holes		6	ACSIEGHY21000
Hy21° bedpan seat bottomless		6	ACSIEGHY21SF0
Hy21° bedpan seat full bottom		6	ACSIEGCUVE000
Seat round, open front, less cover		1	CHSIROOUSC001
Pail with holes		1	CHCUVEHYAL000
Pail with NO holes		1	CHCUVEHYAL001
REPLACEMENT PARTS	PRICE	UNITS PER CASE	PRODUCT CODE
Right arm bracket. Hardware included.		1	PIBLOCBARRDRO
Left arm bracket. Hardware included.		1	PIBLOCBARRGAU
¼"-20 ½" long self-tapping screw for arm bracket. Order 2 units per bracket.		1	PIVISAUTOL205
Height adapter. Hardware included.		1	PI-ADAPHAUT000
Swivel caster. Hardware included.		4	PIROULSA000
Screws kit for caster containing 2 bolts and 2 washers. Order 1 kit per caster.		1 kit	PIKITSROUE8X12N

Parts should be ordered through Hygie at po@hygie.com.

Minimum order for replacements parts: € 30.00

Prices are subject to change without notice.

For questions info@hygie.com